



CROYDE MEDICAL

Maternity Equipment Specialists

CROYDE EXAMINATION COUCHES

The Croyde Examination Couches are built perfectly for general examination and assessment.



The Croyde Examination Couch is available in 2 configurations; 2 Section & 3 Section.

Both models have an excellent low height position to facilitate easy entry & exit, whilst also retaining a high height position that allows it to be set-up for taller clinicians. This feature also provides easy access for cleaning.



Available with a choice of electric or hydraulic profiling functions and a range of upholstery colours, the couches can be ordered specifically to suit the end user.



2 SECTION COUCH

The Croyde 2 Section Examination Couch is strong, stable and smooth in operation during examination and assessment.



3 SECTION COUCH

The Croyde 3 Section Examination Couch allows the patient to be sat up in a chair-like position when longer periods of occupancy are required.



RETRACTABLE WHEEL BASE

The couch is secured into position by activating the brake and lowering the feet onto the floor using the retractable wheel base system.

When the brakes are removed, small castors are deployed to allow the couch to be moved for cleaning purposes.

If the couch needs to be moved between rooms, mobile castors can be ordered instead of the retractable wheel base. This makes moving the couch much easier over longer distances

[CLICK HERE FOR INFORMATION ON HOW TO BUY THE 2 SECTION EXAMINATION COUCH](#)

Click the links below to access the

[Croyde Knowledge Hub](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour Choices](#)

[CLICK HERE FOR INFORMATION ON HOW TO BUY THE 2 SECTION EXAMINATION COUCH](#)

Click the links below to access the

[Croyde Knowledge Hub](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour Choices](#)



CROYDE MEDICAL

Maternity Equipment Specialists

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit 7, Regal Way, Faringdon, Oxon, SN7 7BX

