

CROYDE AVE 2 BIRTHING BED

The AVE 2 Birthing Bed facilitates a wide range of positioning options to optimise labour and delivery when mothers and birthing people need to be on, or near the bed.



Prioritizing both functionality and comfort, the AVE 2 Birthing Bed helps enhance the birthing experience for all involved, fostering a positive and empowering atmosphere in the delivery room.

Every aspect of the bed has been designed with ease of use in mind whilst offering countless positioning options for any birthing person.

The adjustability of the bed also helps encourage birth partner engagement allowing them to be more actively involved and supportive.

The onboard CPR release is also easily accessible in emergency scenarios.



INTEGRATED FOOT SECTION

The unique foot section design means that the foot section can be lowered and stowed underneath the bed, removing the need to lift the foot section off, remove the mattress, or disconnect any parts.



LITHOTOMY LEG RESTS

The integrated leg rests are easily adjusted to find the best position for the mother. The leg rests are folded over and the calf rests can be set into position. There's no need to remove or reattach to the bed.




INDEPENDENT PELVIC TILT

The pelvic tilt allows the seat section to be tilted upwards to prevent the woman slipping down the bed. This can be easily positioned and can also be used to facilitate suturing and vaginal examination.

Upright and active labour is associated with improved outcomes through optimised fetal positioning and movement through the pelvis.
Optimising labour & delivery positions is just as important when mothers and birthing people need to be on, or near the bed.
A birthing bed is designed to help women achieve their most comfortable positions. Twelve popular examples are demonstrated below.

Other equipment may be used to achieve these positions. Consult your midwife.



POSITIONING OPTIONS

The AVE 2 Birthing Bed can be used to optimise positioning during labour & delivery.

Birthing aids such as a beanbag or peanut ball can be used to support positioning.

We partnered with All4Maternity to create a poster showcasing some of the positioning options available.

[DOWNLOAD YOUR POSTER](#)

EVERYTHING YOU NEED TO KNOW

[CLICK HERE FOR INFORMATION ON
HOW TO BUY](#)

Click the links below to access the Croyde Knowledge Hub:

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Videos](#)



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

