



Birth slings give mothers flexible options for vertical hanging and pulling during labour.

With our ceiling-mounted hooks, you can attach slings around the birth or delivery room, and support her in finding comfortable positions.



CEILING SLIDE BAR

The Ceiling Slide Bar is a sliding unit where you can attach a birthing sling at different points. The fabric sling is securely knotted onto a steel ring. The ring attaches to the slide bar via a robust and safe snap hook.

You can install ceiling slide bars above birthing beds and pools. This helps you maximise active birth positioning opportunities for mothers.



STATIC CEILING HOOK

The Static Ceiling Hook has a suspension bar for attaching birth ropes to the ceiling.

Fabric birthing slings are securely knotted onto a steel ring. This attaches to the suspension bar via a robust and safe snap hook system. You have a choice of length and colour of the sling to ensure it suits the design of your delivery room.



EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour Options](#)

[CLICK HERE FOR INFORMATION ON HOW TO BUY](#)



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

