

CROYDE BIRTHING COUCH

Support natural birth with the Croyde Birthing Couch – a birthing couch, seat and overnight double bed, all in one.



The Croyde Birthing Couch helps women find their most comfortable positions during labour and delivery. It supports natural and active childbirth with three setups; unfolded couch, folded over seat or base or queen sized double bed.

The couch is easy to move around the room as it comes with a transport base and handles.

One person can set it up alone. It's ergonomically friendly, with minimal lifting required.

FLEXIBLE DESIGN

The Croyde Birthing Couch helps women find their most comfortable positions during labour and delivery. It supports natural and active childbirth with three setups:



UNFOLDED COUCH

The unfolded, full-length birthing couch provides additional space so the mother can find her most comfortable position. It also encourages engagement with the birth partner.




FOLDED COUCH

The folded couch delivers a comfortable, supportive seat or kneeling base, with additional positions for the midwife or partner. It provides standing labour support through extra height.



DOUBLE BED

The couch can transform into a queen-size double bed for spacious postnatal resting or sleeping. To set up, simply insert the positioning wedge into the base of the couch.



INFECTION CONTROL

The birth couch is covered with two high-quality, medical-grade fabrics. One internal cover protects the foam, and the outer cover encases them both. They are waterproof, breathable and have two-way stretch.

All zips are water repellent and protected by double waterfall flaps that repel fluids. The zips and light underside of the covers make it easy to inspect internal areas.

PU taping on the seams resist fluid ingress, even under intense strain. The taping is monolithically bonded into the material during the manufacturing process making it easy to clean.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the
Croyde Knowledge Hub
[How to Buy](#)

[Dimensions & Specifications \(CBC\)](#)

[Dimensions & Specifications \(CBC+\)](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Videos](#)

SPECIFICATION OPTIONS

CROYDE BIRTHING COUCH - Birthing couch with folding capability, transport base with handles, kneeling mat and positioning wedge. The optional storage cover can have any logo or text printed onto it.

CROYDE BIRTHING COUCH PLUS - Includes the same specification as above but adds a Croyde Birthing Beanbag and a Croyde Positioning Pillow.



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

