

CROYDE BREAST MILK FRIDGES

A selection of milk fridges in a range of sizes to suit your needs.



The Croyde Breast Milk Fridges are available in a range of sizes with 30 to 400 litres of capacity available.

Smaller cabinets are designed to fit underneath or on top of work surfaces. Larger cabinets are equipped with rear wheels for ease of movement.

All of the fridges require a key to unlock. Keeping breast milk secure can contribute to a hospital's Unicef Baby Friendly Initiative accreditation.



REMOTE ALARM SYSTEM

As an optional extra, Croyde Breast Milk Fridges can be wired to a remote alarm (buzzer, auto-dialler, etc) if the set temperature deviates or the power shuts off.



EASY CLEANING

Each fridge is fitted with a painted aluminium interior for easy cleaning. The painted aluminium interior and painted steel exterior are compatible with hospital approved wipes and chlorine-based cleaning solutions.



INTEGRATED DOOR ALARM

If the fridge door is left open or the set temperature deviates, audible and visual alarms signal that the temperature is compromised and the safety of the breast milk is at risk.



PRECISE TEMPERATURE CONTROL

All Croyde Breast Milk Fridges are equipped with digital temperature controls.

They have a temperature range of 1-5°C with a single press 3°C setting for the optimum breast milk temperature.

The quiet internal fan optimises cool air distribution around the chamber.

The painted aluminium interior delivers consistent temperature regulation to keep breast milk cool when the fridge door is open.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Warranty & Servicing](#)

[Cleaning Instructions](#)



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA