

CROYDE CO-SLEEPER COT

Position baby safely next to the bed with the Croyde Co-Sleeper Cot, our hospital baby crib with a home-from-home appearance.



The Croyde Co-Sleeper Cot facilitates closeness and provides instant access to baby from the side of the bed.

This safe and secure cot has height adjustable side rails and lipped edges.

With a light wooden finish and homely design, the Croyde Co-Sleeper Cot has a soft and non-clinical appearance.

Available in two widths with the larger cot being suitable for twins.

ADJUSTABLE SIDE RAILS

The height adjustable side rail protects the baby from rolling out of the cot.

Simply pull up the rail and snap it in the locked position during transport or when parked.

Even with the side rail lowered, the baby is protected from rolling out by the lipped edge.





MATTRESS & SIDE RAIL BUMPER PADS

The Croyde Co-Sleeper Cot is supplied with a washable mattress and separate bumpers.

The bumper pads are easily detached if preferred and the mattress is easily removed for cleaning.

WHEELS & PROTECTIVE BUMPERS

The large, robust rubber rollers have non-tilting arms for easy and almost silent pushing.

The roll extensions serve as stable ram protection to protect the bed and baby in the event of a collision.

The wide base also provides exceptional stability.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

DIMENSIONS

Standard cot dimensions:

Width: 980mm

Depth: 520mm

Height: 915mm

Bed surface: 810 x 380cm

Weight: 14.5kg

Twin cot dimensions:

Width: 980mm

Depth: 680mm

Height: 915mm

Bed surface: 810 x 530cm

Weight: 16kg



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

