

# CROYDE DIASYS EXAMINATION LIGHT

High-quality illumination for examination and suturing.



The Croyde Examination Light is a high-quality light delivering excellent illumination during complex suturing and examination.

The Motion-Matic sensor allows the adjustable LED light intensity to be changed with a hands-free gesture and automatically turns off after 30 mins of inactivity to save energy.

With a stable base and sturdy spring arm, the light head stays perfectly in position when set. The upright has an integrated push handle and large castors ensure movement between rooms is both easy and balanced.



## HIGH-QUALITY ILLUMINATION

The light head includes 4 \* LED lamps delivering excellent illumination of up to 65,000 LUX with a colour rendering index of 96.

This enhanced level of light quality is capable of supporting more complex suturing tasks to avoid the need for patient transfers, specifically to Labour Ward or Theatre.

The robust spring arm design means that it retains the light head position without fear of drifting.



## MOTION-MATIC SENSOR TECHNOLOGY

The innovative sensor technology allows the brightness to be adjusted by simply waving your hand past the sensor.

The intelligent sensor automatically switches the light on when it is pulled into position and turns it off when it is turned upwards. There is no on/off switch.

The sensor also recognises periods of inactivity and automatically turns off after 30 mins. This saves energy and improves the life span.



### MOVING THE LIGHT

Easy-roll castors and a lightweight base make for easy transportation. It is well balanced due to the base design and also has an integrated handle to make pushing more ergonomic. The mains lead can also be wrapped around the handle for safe storage during transit.



### INFECTION CONTROL

The light head is smooth, rounded and easy to clean with a solid back. There are minimal joins in the casing and no on/off switch or button. The whole light is compatible to be cleaned with hospital wipes and chlorine-based solutions.

### EVERYTHING YOU NEED TO KNOW

Click the links below to access the  
Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

### [CLICK HERE FOR INFORMATION ON HOW TO BUY](#)



CROYDE MEDICAL

## CONTACT US

Tel: 01367 241194

Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)

[\*\*www.croydemedical.co.uk\*\*](http://www.croydemedical.co.uk)

Croyde Medical, Unit A4 GBP1, Grove Business Park,  
Downsview Road, Wantage, OX12 9FA

