

# CROYDE UNIVERSAL BIRTH POOL STEPS

The unique design of the Croyde Universal Birth Pool Steps makes them safe and easy to use with any birth pool.



The Croyde Universal Birth Pool Steps have been thoughtfully created with both the birthing person and hospital staff in mind to create a safe and fit for purpose solution that can be used with any birthing pool.

Designed to be used safely with any birth pool, the double step design provides easy access to the pool and the handrail gives additional support & reassurance.

Made from Linear Low Density Polyethylene (LLDPE), the steps are exceptionally robust and highly resistant to damage.



## EASE OF MOVEMENT & STORAGE

The folding handle means that the steps can be conveniently folded to the size of a small suitcase and stored easily in the birth room.

This design means the steps can be either pushed or pulled into position. They can also be tipped back onto their wheels to be moved around.



## STABILITY & BALANCE

The Croyde Stability Bar has been added underneath the base as a counter-balance to ensure stability.

This extra weight lowers the centre of gravity, making it secure and tip-resistant. The base is also flared to maximise the surface area.



### VERSATILE & CLOSE

The uniquely shaped handrail means that the steps can be placed closer to the edge of the pool to reduce the amount of stretching required to reach the rail.

Thoughtfully chosen, the smaller circumference of the rail means that it is comfortable to hold, even for small or swollen hands.



### CLEANABLE & SUSTAINABLE

The steps are fully waterproof and resistant to corrosion and fluids. They can be cleaned with chlorine based solutions and hospital-approved disinfectant wipes.

They are made from 100% recyclable materials, reflecting our commitment to sustainability.

### EVERYTHING YOU NEED TO KNOW

Click the links below to access the  
Croyde Knowledge Hub

[How to Buy](#)

[Dimensions & Specifications](#)


[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Videos](#)

**Weight - 13.5kg**



CROYDE MEDICAL

## CONTACT US

Tel: 01367 241194

Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)

[www.croydemedical.co.uk](http://www.croydemedical.co.uk)

Croyde Medical, Unit A4 GBP1, Grove Business Park,  
Downsview Road, Wantage, OX12 9FA

