

# ELIOS PHOTOTHERAPY WITH INTEGRATED RADIANT WARMING

A phototherapy lamp with integrated radiating warming, featuring an additional built-in examination light.




The Elios delivers high-intensity light that meets the requirements for double phototherapy, ensuring faster and more effective treatment. Its Blue Power LEDs can also be used alongside the white halogen lamps, allowing for clinical examinations to take place during phototherapy sessions.

These light rays can penetrate Perspex, making it ideal for use in closed cots and cribs without interrupting therapy. Ultimately, by providing consistent, uninterrupted phototherapy, the Elios can help to reduce hospital stays and improve patient turnover, enhancing care efficiency.



## HEATING

The 350w heating element delivers radiant warmth from above a crib, changing table or baby bed.



## EXAMINATION

The Elios is equipped with four halogen lamps for excellent patient surface illumination and examination.



## PHOTOTHERAPY

The Elios' six Blue Power LEDs provide consistent phototherapy treatment.

## TOUCHSCREEN CONTROL

The Elios is controlled by its large touchscreen display.

With controls built into the lamp, settings changes happen easily and efficiently

Treatment time is displayed on the screen to easily keep track of a patient's time under the phototherapy lamp between assessments.





## CONSIDERED DESIGN

The Elios has a 1500mm gap between its base and lamp, adjustable to 1800mm.

This clearance means it can sit above any incubator, cot or changing table at any height.

It also has three free moving castors, each with a brake, for ease of movement and stability when stationary.

The wide base gives the Elios extra stability, whilst its shape allows it to slide under cots or beds with ease.

## EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub


[How to Buy](#)

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Warranty & Servicing](#)

[Cleaning Instructions](#)



CROYDE MEDICAL

## CONTACT US

Tel: 01367 241194

Email: [info@croydemedical.co.uk](mailto:info@croydemedical.co.uk)

[\*\*www.croydemedical.co.uk\*\*](http://www.croydemedical.co.uk)

Croyde Medical, Unit A4 GBP1, Grove Business Park,  
Downsview Road, Wantage, OX12 9FA