

VARIOTRAC ROPE SYSTEM

Give women flexibility with the Variotrac Rope System, a dynamic overhead rope sling.



Variotrac provides you with several options for positioning a birthing sling. It enables more movement and convenience within your birth room.

The swivel system can be mounted to a ceiling, wall or floor with 5 possible hook positions.



360° SWIVEL ARM

The Variotrac system can rotate through 360° and lock securely at 15° increments.

You can position it at different points in your delivery room. This gives women something to pull down on and also can be used to support women when entering or exiting a birth pool.



COMPATIBLE WITH BIRTHING EQUIPMENT

You can use Variotrac with other birth aids to make childbirth more manageable. For instance, it works well with our birthing balls, birth stools, beds, couches and pools.

The Variotrac Rope System functions perfectly with the Duoetto birthing system, Pelviball and Dullstein birth stool



MOUNTING OPTIONS & TECHNICAL SUPPORT

We can install the Variotrac Rope System into any birth or delivery room. It suits all spaces as we can mount the system to your ceiling, wall or floor.

The ceiling-mounted version can fit directly onto the ceiling. Alternatively, it can attach via an intermediate construction for suspended ceilings. The centre of the mount to the sling's eyelet is 1,300mm, but we can provide a 1,700mm arm if needed.

Our technical team is happy to provide you with CAD drawings bespoke to your room.

EVERYTHING YOU NEED TO KNOW

Click the links below to access the Croyde Knowledge Hub

[Dimensions & Specifications](#)

[Delivery & Implementation](#)

[Cleaning Instructions](#)

[Warranty & Servicing](#)

[Colour Options](#)

[CLICK HERE FOR INFORMATION ON HOW TO BUY](#)



CROYDE MEDICAL

CONTACT US

Tel: 01367 241194

Email: info@croydemedical.co.uk

www.croydemedical.co.uk

Croyde Medical, Unit A4 GBP1, Grove Business Park,
Downsview Road, Wantage, OX12 9FA

